



LUSCIIAPP
Peripheral Devices:
Recording blood
oxygen levels



Luscii
an OMRON Healthcare Service

This guide is supplementary to the manufacturer(s) Instructions for Use (IFU) and does not replace the IFU issued by them. Graphnet is acting as a distributor only. Please refer to the manufacturer's IFU which should accompany the device. For the Luscii App please refer to the IFU available on their website. www.luscii.co.uk

Blood Oxygen Monitor Devices

Your finger clip SpO2 / blood oxygen monitor will look like one of these (or similar), it is used to measure the oxygen in your blood.

Taking blood oxygen levels

1. Rest for at least 15 minutes before taking blood oxygen levels (SpO2)
2. When the question arises that prompts for blood oxygen levels, pick up the blood oxygen monitor that has been provided
3. Place the blood oxygen monitor on the middle finger of either hand as instructed. The monitor will automatically display the measurement






Recording blood oxygen levels

4. To accurately enter the measurements:

Select how you would like to enter your results, click next.

Scroll up or down to find your oxygen level as displayed on the blood oxygen monitoring device, click next.

When you are satisfied that this is your reading, click done (iPhone) or Next (Android) to submit.





Graphnet 
Transforming Care

 **Luscii**
an OMRON Healthcare Service

Graphnet Health Limited
Unit 4A/B Newton Court
Farraday Road
Wavertree Technology Park
L13 1EJ

T: +44 (0)3330 771 988
E: liverpoolsupport@graphnethealth.com
www.graphnethealth.com