



Luscii
an OMRON Healthcare Service


LUSCIIAPP
Peripheral Devices:
Recording temperature

This guide is supplementary to the manufacturer(s) Instructions for Use (IFU) and does not replace the IFU issued by them. Graphnet is acting as a distributor only. Please refer to the manufacturer's IFU which should accompany the device. For the Luscii App please refer to the IFU available on their website. www.luscii.co.uk

Recording temperature



You will have been provided with an infrared thermometer.

- When using the thermometer for the first time, move the battery's insulating piece away.
- Press the power button to power on the thermometer.
- When measuring adults, place the switching device in the "adult" mode. When measuring children, place the switching device in "child" mode and the corresponding symbol is displayed on the screen.





Using the correct modes on your thermometer

- When taking your temperature, the display should look like this with a symbol of an adult head and the measurement in °C
- If you see either of these symbols, press the MODE button on the side of the device until you see the adult head.
- If the display is showing °F you need to change to °C



To do this, with the device switched on, press and hold the measurement button for 3 seconds until the display changes.



Infrared / Forehead Thermometer

- Press the on/off button to activate the thermometer
- The thermometer is ready to use when an indication icon appears on the display. This will vary depending on the brand and type of thermometer.
- Point the infrared end directly towards your forehead, between the brows and approximately 5cm / 2 Inches away.
- Press the start button which may produce a beep sound. The measurement will be displayed on the screen.
- Write the measurement down as the thermometer will power off automatically after around 10 seconds.



× **Temperature**

Enter your temperature

35	
36	.0
37	.1

°C

Back
Next

Recording Temperature on LUSCIIAPP

- Scroll up and down the numbers on the left until you reach your temperature from the thermometer and then do the same for the decimal numbers on the right.
- When you have reached your reading, click next
- If the reading on this screen matches your reading from the thermometer, click Done (iPhone) or Next (Android) to submit your recording



Graphnet 
Transforming Care



Luscii
an OMRON Healthcare Service

Graphnet Health Limited
Unit 4A/B Newton Court
Farraday Road
Wavertree Technology Park
L13 1EJ

T: +44 (0)3330 771 988
E: liverpoolsupport@graphnethealth.com
www.graphnethealth.com