



Luscii
an OMRON Healthcare Service


LUSCIIAPP
Recording your ECG

This guide is supplementary to the manufacturer(s) Instructions for Use (IFU) and does not replace the IFU issued by them. Graphnet is acting as a distributor only. Please refer to the manufacturer's IFU which should accompany the device. For the Luscii App please refer to the IFU available on their website. www.luscii.co.uk

Recording your ECG

In your list of actions, click on 'ECG'

Click on 'Single-Lead ECG'





Ensure the AliveCor device is within 30cm from your smart device

Rest two or more fingers on the AliveCor device

Rest your arms on a flat surface and remain still with your fingers on the electrodes for 30 seconds

After 30 seconds, your ECG reading will be complete. Click the green tick.



Your ECG results are now ready. Click on Done.



For technical assistance call 0151 294 3364

Graphnet 
Transforming Care

 **Luscii**
an OMRON Healthcare Service

Graphnet Health Limited
Unit 4A/B Newton Court
Farraday Road
Wavertree Technology Park
L13 1EJ

T: +44 (0)3330 771 988
E: liverpoolsupport@graphnethealth.com
www.graphnethealth.com